



# Hardness Test Blocks

<b>Knoop / Vickers Test Blocks</b>	
<b>Catalog Number</b>	<b>Description</b>
9200-200	Test Block, Microhardness, Knoop 10 to 499 grams- 1-25" dia solid metal - In accordance with ASTM E-384 <i>*Specify gram load and nominal hardness value at time of order</i>
9200-201	Test Block, Microhardness, Knoop 500 to 1000 grams- 1-25" dia, solid metal - In accordance with ASTM E-384 <i>*Specify gram load and nominal hardness value at time of order</i>
9200-202	Test Block, Microhardness, Vickers 50 to 499 grams- 1-25" dia, solid metal - In accordance with ASTM E-384 <i>*Specify gram load and nominal hardness value at time of order</i>
9200-203	Test Block, Microhardness, Vickers 500 to 1000 grams- 1-25" dia, solid metal - In accordance with ASTM E-384 <i>*Specify gram load and nominal hardness value at time of order</i>
900020515	Test Block, MacroVickers 1kg to 1000kg - In accordance with ASTM E-92 <i>*Specify gram load and nominal hardness value at time of order</i>
900020511	Additional Calibration on Microhardness Test Block-Specify nominal hardness value and load when ordering
9200-100	Vickers Microhardness Test Block, between 500 and 800 HV, 500 grams- Solid, non-mounted block, 1-25" dia-
9200-104	Knoop Microhardness Test Block, between 500 and 800 HK, 500 grams- Solid, non-mounted block, 1-25" dia-
9200-105	Vickers, 765 HV 10 kg, Hardness Test Block

Americas:  
Tel: 1-800-695-4273  
hardness\_inquiry@instron.com

Europe:  
Tel: 0049 - 241 997 2800  
info@wolpertwilson.com

Asia:  
Tel: 0086 - 21 6476 6110  
Asia.sales@wolpertwilson.com